

cat. #1021-500F



500 reactions

## EPPiC Fast

Enzymatic Post-PCR immediate Cleanup



A&A BIOTECHNOLOGY, Aleja Zwycięstwa 96/98, 81-451 Gdynia, Poland  
tel. +48 58 6982194, fax +48 58 6228578, www.aabiot.com, info@aabiot.com

## EPPiC Fast

Enzymatic Post-PCR immediate Cleanup

### Quick and effective PCR fragment purification

The EPPiC Fast mixture contains two enzymes that effectively degrade dNTPs and primer left-overs from previous PCR mixtures while leaving the double stranded DNA PCR products untouched.

The EPPiC Fast enzymes are active at 37 °C in standard buffers used in PCR and are completely thermally inactivated by 1 min at 80 °C incubation. The EPPiC Fast purified PCR fragments are free from dNTPs and original PCR primers and may be directly used in downstream applications including nested, second round PCR, SNP analysis, molecular cloning and cycle sequencing reactions.

#### The EPPiC Fast mixture contents:

- recombinant thermolabile nucleotide hydrolase
- recombinant exonuclease I with increased efficiency

#### Store at -20 °C

*Note: Unlike the other enzymatic mixtures used for PCR fragment clean up, EPPiC Fast mixture does NOT remove the 5'-phosphate groups from PCR products obtained with phosphorylated PCR primers. Therefore subsequent cloning of EPPiC-purified, phosphorylated PCR products does NOT require extra phosphorylation of 5'-ends.*

500 reactions



cat. # 1021-500F

## EPPiC Fast

Enzymatic Post-PCR immediate Cleanup

### Purification protocol of PCR product:

#### To prepare the samples we recommend using a thermocycler.

1. Briefly spin the EPPiC Fast mixture and place the tube on ice.
2. Add 2 µl of the EPPiC Fast mixture per each 10 µl of post-PCR mixture. Mix by pipetting.  
*Note! For example use 1.0 µl of the EPPiC Fast per 5 µl of PCR or proportionally 4 µl of the EPPiC Fast per 20 µl of PCR.*
3. Place the tube in the thermocycler and start the programme: 37 °C - 5 min, then 80 °C - 1 min to inactivate the EPPiC Fast mixture.
4. Briefly spin the sample.
5. The resulting PCR product is ready for downstream application and may be stored at -20 °C until needed.

For R&D use only

500 reactions



cat. # 1021-500F

## EPPiC Fast

Enzymatic Post-PCR immediate Cleanup

### Related products:

<b>EPPiC:</b>	
100 reactions	cat. # 1021-100
500 reactions	cat. # 1021-500
<b>Clean-Up:</b>	
50 reactions	cat. # 021-50
250 reactions	cat. # 021-250
<b>Clean-Up Concentrator:</b>	
50 reactions	cat. # 021-50C
250 reactions	cat. # 021-250C
<b>Gel-Out:</b>	
50 reactions	cat. # 023-50
250 reactions	cat. # 023-250
<b>Gel-Out Concentrator:</b>	
50 reactions	cat. # 023-50C
250 reactions	cat. # 023-250C
<b>Exonuclease I:</b>	
1000 U	cat. # 1020-1
5000 U	cat. # 1020-5

### PCR Master Mixes (200 reactions in 25 µl):

PCR Mix	cat. # 2005-100
PCR Mix Plus	cat. # 2005-100P
2CR Mix Plus HGC	cat. # 2005-100G
StartWarm HS-PCR Mix	cat. # 2017-100
PCR Mix Rapid	cat. # 2009-100

500 reactions



cat. # 1021-500F